

AUTOCLAVE STEAM STERILIZER

- Fully Automatic Autoclave GMP Compliant Series Lab Type (Available) Semi Automatic, Standard



- Full Automatic Laboratory Autoclave Radial Locking System & Paddle Free System



Application

Micro Biology Labs, Pharma / Biotech R & D Centres & Drug Discovery etc.

Options

IQ/OQ/PQ Documentation. CCR - Circular Chart Recorder. DPQ - Digitak Pressure Display. Di4-4 Channel Data Logger. SCR - Strip Chart Recorder.

Features

- Unique Paddle Operated Design
- Lower Working Height
- Microprocessor Based Precise Controlling
- User Suitable Multi Program Cycle
- Automatic Purging, Venting

SPECIFICATIONS :

Model	AC-35-FA	AC-52-FA	AC-111-FA	AC-178-FA
Capacity	35 ltrs	52 ltrs	111 ltrs	178 ltrs
Chamber Dimension	300 x 500 mm	350 x 550 mm	450 x 700mm	550 x 750 mm
Construction	SS 304	SS 304	SS 304	SS 304
Power Consumption	3 Kw	3 Kw	4 Kw	6 Kw
Operating Pressure	15 PSI	15 PSI	15 PSI	15 PSI
Operating Temperature	121°C	121°C	121°C	121°C
Sterilization Time	0-99 mins.	0-99 mins.	0-99 mins.	0-99 mins.
Voltage	230 V	230 V	230 V	230 V
Control	Digital	Digital	Digital	Digital
Shipping Size (mm)	780x570x1035	785x630x1055	1050x875x1425	1125x975x1475