

PORTABLE AUTOCLAVE SEMI AUTOMATIC

SPECIFICATIONS :



- Operating System: semi-automatic Table Top Autoclave
- MOC : Complete Stainless Steel
- Inner body heavy gauge SS304
- Outer shell SS304
- Temperature Control: Temperature is controlled with the help of Digital PID Controller in conjunction with PT-100 RTD Temperature sensor;
- Fitted with pressure gauge, safety valve and steam release valve
- Electrically operated on 230V AC/ 50Hz.

ORDERING INFORMATION :

Cat. No.	Size of Working Chamber (inch)	Size of Working Chamber (cm)	Volume Ltrs	Rating (KW)
AC - 21 - SA	12 x ht. 12	30 dia x 30 ht	21	2.0