

BUNSEN BURNER - DIGITAL



A **DIGITAL BUNSEN BURNER** is used for non-contact heating of test tubes, crucibles, flasks, and beakers. It can also be used to sterilize inoculating needles for microbiological work.

TECHNICAL SPECIFICATIONS :

- Conically shaped refractory crucible
- **Size of cup** : 55 mm Dia x 47 mm Ht.
- Stainless Steel cylindrical body SS 304 with mirror finish
- Control Panel made of SS 304 with mirror finish
- Maximum working temp. 450 °C. Approximate Temperature.
- The temperature is controlled by digital PID temperature controller.
- Fuse fitted on backside for overload protection. (SAFETY)
- **Power** : 350 W
- Works at 1 Phase 230 Volts 50 Hz

CAT. NO.	Model	Type	Shape	Watts
BBD - 001	BUNSEN BURNER - Digital	Refractory Crucible	Conical	350

HSN CODE – 8419 1920 GST % - 18 %