

## HEATING MANTEL ANALOG / DIGITAL



**HEATING MANTLES** are used for heating or tempering organic liquids placed in reaction kettles, round-bottomed flasks, or relevant reaction vessels required for the boiling, evaporation, distillation, or extraction process.

### HEATING MANTLE

- MS Powder Coated Aluminium metal Body.
- Energy Regulator, Switch & Indicating Lamp
- **Temperature range** : (Ambient + 5) to 350 ° C
- Max. Temperature upto 350 °C Aprox.
- Fuse fitted on backside for overload protection. (SAFETY)

### ORDERING INFORMATION :

| CAT. NO.               | MODEL                           | CAPACITY (ML) | WATTS   |
|------------------------|---------------------------------|---------------|---------|
| HM - 0002              | Heating Mantel Analog           | 100           | 60      |
| HM - 0003              | Heating Mantel Analog           | 250           | 85      |
| HM - 0004              | Heating Mantel Analog           | 500           | 150     |
| HM - 0005 / HMD - 0005 | Heating Mantel Analog / Digital | 1000          | 200     |
| HM - 0006 / HMD - 0006 | Heating Mantel Analog / Digital | 2000          | 300     |
| HM - 0007 / HMD - 0007 | Heating Mantel Analog / Digital | 3000          | 450     |
| HM - 0008 / HMD - 0008 | Heating Mantel Analog / Digital | 5000          | 2 X 300 |
| HM - 0009              | Heating Mantel Analog           | 10000         | 2 X 600 |
| HM - 0010              | Heating Mantel Analog           | 20000         | 3 X 665 |

HSN CODE – 8419 1920 GST % - 18 %