

BENCH TOP ORBITAL SHAKER - HEAVY DUTY



SPECIFICATIONS :

- **Body Material** - MS Powder coated
- **Tray Material** - Stainless steel SS 304
- **Shaking Speed** : 20 - 250 Orbital strokes per min (without Load)
- Works on 230 Volt AC Single phase supply.

SPECIAL FEATURES :

- Brushless Induction Motor, for 24 hrs operation.
- DC DRIVE
- **Digital Voltmeter** : With LCD display for better view in the dark room.

ORDERING INFORMATION :

CAT. NO.	Tray Size in mm	RPM SPEED	Digital Speed Indication
OSR - 81 100	435 x 435	20 - 250	Yes
OSR - 81 750	750 x 750	20 - 250	Yes

Note : SS Lotus clamps for holding the flasks in the shakers to be order separately.

CAT. NO.	Flask Capacity (ml)	Platform size 750 mm x 750 mm can accommodate the following combinations	Platform size 435mm x 435 mm can accommodate the following combinations
SS CLAMP - 501	50	49	25
SS CLAMP - 502	100	49	25
SS CLAMP - 503	250	49	16
SS CLAMP - 504	500	25	16
SS CLAMP - 505	1000	25	8
SS CLAMP - 506	2000	9	4

HSN CODE – 8479 8200 GST % - 18 %