

DIGITAL VDRL SHAKER (Premium Model)

Microprocessor Based



FEATURES :

- Adjustable speed control provides gentle to vigorous agitation.
- Microprocessor control with simultaneous display of time and speed
- Brushless DC motor ensures maintenance-free, reliable and quiet operation Audio alarm signals program completion

SPECIFICATIONS :

• Model No.	DVS - 304
• Mounting Type	BenchTop
• Controller	Microprocessor
• Type Of Motor	Brushless
• Stepless Frequency Speed Control	Yes
• Shaking Speed Range	1-180 RPM (1 to 300 RPM Optional)
• Shaking Amplitude	20mm
• Compact Counter Balance Drive Mechanism	Yes
• Shaking Motion	Orbital
• Platform Size in mm	304mm x 304mm
• Display Type	LCD
• Auto restart after power failure	Yes
• Timer Range	999h 59 min digital timer
• Operation	Continuous / Timer based
• Display	16 x 2 Numeric LCD
• Max Load Capacity	2 Kg
• Power Consumption	20 W
• Power Supply	220 VAC- 50Hz
• Weight	10 Kg

HSN CODE - 8479 8200 GST % - 18 %