

ORBITAL SHAKER (Premium Model)



FEATURES :

- Adjustable speed control provides gentle to vigorous agitation
- Microprocessor control with simultaneous display of time and speed
- Brushless DC motor ensures maintenance-free, reliable and quiet operation
- Audio alarm signals program completion
- ABS (non broken type food grade) make lotus clamps

SPECIFICATIONS :

- Shaking Motion: Orbital
- Shaking orbit: 20 mm
- Shaking Speed Range: 01 - 300 RPM
- Timer Range: 99 h 59 min digital timer
- Operation: Continuous / Timer based
- Display: 16 x 2 Numeric LCD
- Overall dimension: W - 410 x L - 530 x H - 200 with Lotus Clamp
- Power Consumption: 10 W
- Power Supply: 220 VAC - 50 Hz

CAT. NO.	Tray Size in mm	Digital Speed Indication	Digital Timer
DOS - 304	304 x 304	Yes	Yes
DOS - 457	457 x 457	Yes	Yes
DOS - 635	635 x 635	Yes	Yes

Note : SS Lotus clamps for holding the flasks in the shakers to be order separately.

CAT. NO.	Flask Capacity (ml)	Tray size 305 X 305 mm	Tray size 456 X 456 mm	Tray size 635 X 635 mm
CLAMP - 501	50	25	25	36
CLAMP - 502	100	25	25	36
CLAMP - 503	250	9	16	25
CLAMP - 504	500	9	16	25
CLAMP - 505	1000	5	9	16

HSN CODE - 8479 8200 GST % - 18 %