

UNIVERSAL CELL MIXER



A **UNIVERSAL CELL MIXER** is a device that can be used for a variety of applications, including:

- **Tissue culture:** Suitable for mini tissue culture work
- **Blood suspensions:** Can be used for blood suspensions without sedimentation or foaming
- **Vigorous agitation:** Can be used for vigorous agitation and end-to-end mixing
- **Mounting tubes:** Tubes can be mounted at any angle from horizontal to vertical

SPECIFICATIONS :

- Suitable for mini tissue culture work and blood suspensions without sedimentation or foaming.
- Also suitable for vigorous agitation and end-to-end mixing.
- Tubes can be mounted at any angle from horizontal to vertical, tubes concerning axis of rotation.
- Different sizes from 12 mm to 25 mm dia can be accommodated
- Speed Variable 10 - 100 RPM (without load).
- Electrically operated on 230V AC single-phase power.

CAT. NO.	Tube Dimensions
CM - 001	UNIVERSAL CELL MIXER with 0.5 ml (Ø 8 mm) Eppendorf tubes - 12 nos & 1.5 ml (Ø 11 mm) Eppendorf tubes - 12 nos
CM - 002	UNIVERSAL CELL MIXER with 15 ml (Ø 13 mm) Eppendorf tubes - 08 nos
CM - 003	UNIVERSAL CELL MIXER with 50 ml (Ø 31 mm) Eppendorf tubes - 08 nos

HSN CODE – 8479 8200 GST % - 18 %