

WATER BATH WITH THERMOSTATIC / DIGITAL CONTROLLER



SPECIFICATIONS :

- Inner chamber is made of thick gauge S.S. 304 grade
- Outer Body with M.S. Powder coating/SS 304
- The gap between the two walls is filled with glass wool.
- Specially made extra long heating elements for high heat dissipation, uniformity & longer life.
- Concentric rings of 3" dia i.e. (7.5 cm) made up of SS 304 with 6 & 12 holes.
- Temperature range: Ambient + 5°C to 99°C
- Temperature controlled by Thermostatic control / Digital PID Controller.
- Provided with indicating lamp & on/off switch
- Electrically operated on 230V AC single phase power

CAT. NO.	Capacity In LTRS	Model	Size In MM	Holes	Inner Chamber	Outer Chamber	Watts
WBD - MS 006	14	STD	325 X 300 X 150	6	SS 304	MS Powder coated	1000
WBD - MS 012	22	STD	500 X 300 X 150	12	SS 304	MS Powder coated	2000
WBD - SS 006	14	GMP	325 X 300 X 150	6	SS 304	SS 304	1000
WBD - SS 012	22	GMP	500 X 300 X 150	12	SS 304	SS 304	2000
WBT - MS 006	14	STD	325 X 300 X 150	6	SS 304	MS Powder coated	1000
WBT - MS 012	22	STD	500 X 300 X 150	12	SS 304	MS Powder coated	2000
WBT - SS 006	14	GMP	325 X 300 X 150	6	SS 304	SS 304	1000
WBT - SS 012	22	GMP	500 X 300 X 150	12	SS 304	SS 304	2000

HSN - 8419 1920 GST % - 18 %